VTYT

To 100 mL VT medium add 10 mg yeast extract and 20 mg tryptone.

Reference

Ichimura, T., Itoh, T. 1977 17. Preservation methods of microalgae (I) [17. Bisaisôrui no hozonhô (I)]. In *Preservation methods of microorganisms [Biseibutsu Hozonhô]*, Ed. by Nei, T., University of Tokyo Press, Tokyo, p. 355-373 (in Japanese without English title).

VT

Ca(NO ₃) ₂ · 4H ₂ O	11.78 mg
β -Na ₂ glycerophosphate · 5H ₂ O	5 mg
MgSO ₄ · 7H ₂ O	4 mg
KCl	5 mg
Vitamin B ₁₂	0.01 μg
Biotin	0.01 µg
Thiamine HCl	1 μg
PIV metals	0.3 mL
Glycylglycine	50 mg
Distilled water	99.7 mL
pH 7.5	

Reference

Provasoli, L., Pintner, I. J. 1959 Artificial media for fresh-water algae: problems and suggestions. In *The Ecology of Algae. Spec. Pub. No. 2,*, Eds. by Tryon, C. A., Jr. & Hartmann, R. T., Pymatuning Laboratory of Field Biology, University of Pittsburgh, Pittsburgh, p. 84-96.

Starr, R. C. 1973 Special methods-dry soil samples. In *Handbook of Phycological Methods. Culture Methods and Growth Measurements*, Ed. by Stein, J. R., Cambridge University Press, Cambridge, p. 159-167.

P IV metals

Na ₂ EDTA · 2H ₂ O	100 mg
FeCl ₃ · 6H ₂ O	19.6 mg
MnCl ₂ · 4H ₂ O	3.6 mg
ZnCl ₂ ¹⁾	1.04 mg
CoCl₂ · 6H₂O	0.4 mg
Na ₂ MoO ₄ · 2H ₂ O	0.25 mg
Distilled water	100 mL

1) In the NIES-Collection, 1.04 mg ZnCl₂ is replaced by 2.2mg ZnSO₄ · 7H₂O.

Reference

Provasoli, L., Pintner, I. J. 1959 Artificial media for fresh-water algae: problems and suggestions. In *The Ecology of Algae. Spec. Pub. No. 2,*, Eds. by Tryon, C. A., Jr. & Hartmann, R. T., Pymatuning Laboratory of Field Biology, University of Pittsburgh, Pittsburgh, p. 84-96.